


The Thermo Scientific Nicolet™  
Almega™ XR family of dispersive  
Raman spectrometers are a series  
of research-grade Raman microscope  
systems designed for laboratory use.  
The highly-configurable, powerful,  
and productive instruments provide  
you with the answers to your  
application problems.

## Nicolet Almega XR

High-Performance Dispersive Raman  
Spectrometers



Dispersive Raman spectrometers utilize a form of vibrational spectroscopy to provide information regarding molecular structure. Raman spectra are used for identification and quantification of both organic and inorganic materials and are very sensitive to minor changes in molecular structure.

### Versatile, Easy Sampling


The Nicolet Almega XR system has many applications in the pharmaceutical, forensic, and semiconductor fields, and features specialized sampling options for many applications in these areas. The system can be configured for both micro and macro sampling of liquids and solids. No sample preparation nor dilution is necessary. Sampling is nondestructive and can often be performed directly through glass or plastic packaging, providing an excellent solution when non-invasive sampling is required. Aqueous-

based solutions are readily analyzed with the Nicolet Almega XR spectrometer. Dispersive Raman provides useful information about these samples, unlike other vibrational spectroscopy techniques. Remote sampling with fiber-optics is also an option.

### Raman Microscopy

The Thermo Scientific Nicolet Almega XR microscope provides exceptional micro-Raman performance. The microscope is interfaced with the main spectrometer entirely using high-throughput, free-space optics. No fiber-optics are used. This provides a very small excitation laser spot size at the sample, eliminating the need to excite excessively large regions of the sample when working with very small samples. This greatly reduces the amount of interfering fluorescence generated by many samples. Combining this feature with the use of the integral confocal

aperture and patented instrument alignment provides exceptional control of fluorescence, superior spatial resolution, and the ability to perform depth-profiling analysis.



Notable features of the Nicolet Almega XR micro-Raman system include:

- Diffraction-limited spatial resolution
- Confocal microscopy and confocal depth profiling
- Polarization studies with automated polarizer option
- Crystal orientation independence with depolarizing kit option
- Selection of precise sample positioning options
  - Manual 2" x 3" stage
  - Motorized 4" x 3" stage
- Well-plate adapter
- Automated alignment to ensure peak performance at all times

### Visual Microscopy

The Nicolet Almega XR system provides excellent visible viewing capabilities in addition to high-performance Raman data collection. A high-quality Olympus BX series visible microscope provides the basis for visual



microscopy on the system. These capabilities are augmented by either a dedicated color-video camera or a high-quality trinocular to provide viewing and image capture on the computer monitor or direct visual viewing for detailed examination. Additionally, the Olympus base permits the Nicolet


Almega XR system to be equipped with a broad range of microscopy viewing options.

Viewing features include:

- Brightfield and darkfield viewing
- Transmission and reflection viewing and illumination option
- Polarized light viewing
- Large selection of standard and long working-distance objectives
- Motorized objective turret on some models
- Autofocus option with motorized stages

### Microscope Mapping

OMNIC™ AtPlus™ software combines the high spatial resolution of the Nicolet Almega XR system with automated data collection to produce extremely high-quality, information-rich Raman area maps. The maps can be displayed as 3D representations or chemical contour plots and each map includes a corresponding video image.




Mapping features include:

- Full spectral information available at every position on map
- Data quality of maps matches that of single point collection
- Completely automated data collection
- Mosaic™ imaging of areas larger than the field of view
- Supports perpetual reanalysis of maps to generate different chemical profiles
- Extensive array of post-processing routines can be applied before generating chemical profiles (e.g., baseline correction, smoothing and spectral subtraction)
- Storage of spectral data, video image, and chemical profiles in one file

### Macro Sampling

Since not all samples are microscopic, the Nicolet Almega XR spectrometer also provides macro sampling capability.




The optional macro compartment includes:

- Entirely software controlled switching between microscope and macro compartment
- Visible and safe laser beam path locator to facilitate optimizing custom accessories
- Automated sampling positioning with the x-y-z stage accessory
- Real-time spectral feedback when optimizing sampling geometries
- Polarization options
- Autosampler option for liquid and solid samples

### Flexible Options

The Nicolet Almega XR system is highly configurable and offers numerous options to meet the needs of specialized applications. Nearly all options can be added to systems in the field at any point. Configure the Nicolet Almega XR system to meet your needs today and have confidence that you will be able to expand your capability in the future.



Flexible options include:

- Wide range of excitation lasers
- Port to bring in external excitation laser
- Simultaneous mounting of two excitation lasers and up to four gratings
- High and low resolution gratings
- Polarizer/depolarizer kits
- Specialized stage options
- Above and below ambient temperature sampling options
- Remote sampling options through fiber optics
- Extensive range of Olympus visual microscopy options

### Automation Makes Your Job Easier

The Nicolet Almega XR system offers an unmatched level of automation. Automation reduces manual operations, which removes subjectivity from the analysis, makes it easy to maintain the instrument at peak performance, and greatly increases productivity.



Automation features and options include:

- Patented† method for system autoalignment
- Automated grating switching
- Automated laser switching and alignment
- Automated Rayleigh rejection filter switching
- Automated calibration and referencing
- Automated white-light correction
- Automated dark current correction for instrument and detector background
- Automated stage control and mapping with motorized stages
- Autofocus option
- Automated switching between microscope and macro compartment

† U.S. Patent 6,661,509



## High-Performance Components

The Nicolet Almega XR system uses extremely high-performance optical components that are selected to be optimal for the spectral range of intended use. High brightness lasers permit focusing to extremely small spot sizes, avoiding excess fluorescence in the area surrounding the region of interest. High brightness lasers also provide energy densities at the sample that are far higher than ordinary lasers, sometimes in excess of 30 times higher, with much lower power outputs. High energy density at the sample enhances the Raman signal, improving results. The performance increases offered by these lasers are especially prominent with very small samples.




High-performance components available include:

- High-brightness lasers
- Gratings optimized for every laser frequency
- High-throughput filters
- High-sensitivity, TE-cooled silicon CCD array detector
- Wavelength optimal reflective optics coatings


## Complete Software Tools

The Nicolet Almega XR system is controlled by OMNIC software. OMNIC is a comprehensive and completely integrated package for data collection, data manipulation, library searching, report generation, and much more. It has been specifically designed to meet the needs of vibrational spectroscopy analysis and has incorporated the feedback of thousands of users to become the standard that it is today. In addition, specialized modules that provide comprehensive toolkits for specific applications are available. All software packages include extensive documentation.



Overview of available software:

- **OMNIC** – for instrument control, data collection, data manipulation, library searching, library creation, and report generation.
- **TQ Analyst™** – comprehensive quantitative analysis and chemometric modeling package supporting Beer's law, CLS, MLR, PCR, and PLS algorithms along with several popular classification algorithms.
- **OMNIC Macros** – provide a simple graphical interface to allow users to capture commonly used software procedures and SOPs as simple push button operations.
- **Atlas** – spectroscopic mapping and chemical profiling software.
- **Array Automation** – automate data collection and analysis for array based samples such as well-plates and blister packs.
- **Data Integrity Software** – data security and digital signatures module for OMNIC.



## Raman Libraries

When trying to identify unknown materials, collecting the spectrum is only half of the job. An extensive collection of high quality reference spectra will greatly enhance your ability to identify unknowns. We provide the largest collection of Raman spectral libraries offered anywhere. All library spectra have been collected by FT-Raman, so they are fluorescence free and cover the full spectral range. The collection consists of over 15,000 spectra, and is still growing. In addition, OMNIC software provides you with the ability to create your own libraries so you can grow your own collection with every new material that you measure.

## Laser Safety

Laser safety is a very important consideration when performing Raman experiments. The Nicolet Almega XR system offers an option to make it a fully integrated Class I laser product.\* Configured with this option the Nicolet Almega XR system can reside in any standard laboratory environment with no need for special laser safety provisions, such as goggles, isolated areas, or interlocked lab entrance doorways.

\* If the Class I option is not selected then the system is a Class IIIb laser product.

## Nicolet Almega XR Dispersive Raman Spectrometer

### System

Olympus microscope, reflection illuminator	Standard
Transmission illuminator	Option
Supports brightfield and darkfield illumination*	Standard
Confocal optics including confocal aperture	Standard
Maintenance-free, extended-range TE-cooled Si CCD	Standard
Automated system calibration	Standard
Fully automated optical alignment	Standard
Class I laser safe enclosure for microscope	Option
Class I laser safe macro compartment	Option
Automated polarization analyzer	Option
Depolarization kit	Option

### Performance

Laser options	785 nm, 780 nm 633 nm, 532 nm, 473 nm (Others available as special order)
Grating options	High and low resolution for each wavelength
Spectral range:	
Absolute	400 – 1050 nm
For 785 nm excitation	100 – 3100 cm <sup>-1</sup> Raman shift
For 780 nm excitation	100 – 3300 cm <sup>-1</sup> Raman shift
For 633 nm excitation	100 – 4000 cm <sup>-1</sup> Raman shift
For 532 nm excitation	100 – 4000 cm <sup>-1</sup> Raman shift
For 473 nm excitation	100 – 4000 cm <sup>-1</sup> Raman shift
Spectral resolution	2 cm <sup>-1</sup> ** FWHM (1 cm <sup>-1</sup> per CCD pixel element)
Spatial resolution	Diffraction limited (1 μm) **
Confocal resolution	2 μm
Confocal aperture	Size selectable slit or pinhole

### Microscope

Stages:	
2" x 3" travel stage	Manual positioning
4" x 3" travel motorized stage	0.1 μm steps, 1.0 μm repeatability
Well plate automation option	Up to 1,536-well plates
Standard working distance objectives	4X, 10X, 20X, 50X, 100X
Long working distance objectives	10X, 20X, 50X, 100X
Microscope sampling geometry	180° Refractive optics

### Physical Dimensions

Optical bench configured with standard microscope	114 cm (W)***, 77 cm (D), 89 cm (H)****
Weight	295 lbs., 134 kg
Optical bench with 300 mm wafer stage	163 cm (W), 89 cm (D) 89 cm (H)**** (Additional clearance needed below front end to open door)
Macro sample compartment dimensions	29 cm (W), 31 cm (D), 28 cm (H)
Macro compartment beam height	8.9 cm

### Warranty

12-month warranty standard on the complete Nicolet Almega XR system. Extended warranties are available.

### Regulatory approval



\* Requires appropriate objective

\*\* Varies with wavelength of excitation laser

\*\*\* Additional 15 cm clearance needed on right side for air cooling vent

\*\*\*\* Some lasers extend beyond the optical bench and will exceed these dimensions

©2007 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Africa** +43 1 333 5034 127  
**Australia** +61 2 8844 9500  
**Austria** +43 1 333 50340  
**Belgium** +32 2 482 30 30  
**Canada** +1 800 530 8447  
**China** +86 10 5850 3588

**Denmark** +45 70 23 62 60  
**Europe-Other** +43 1 333 5034 127  
**France** +33 1 60 92 48 00  
**Germany** +49 6103 408 1014  
**India** +91 22 6742 9434  
**Italy** +39 02 950 591

**Japan** +81 45 453 9100  
**Latin America** +1 608 276 5659  
**Middle East** +43 1 333 5034 127  
**Netherlands** +31 76 587 98 88  
**South Africa** +27 11 570 1840  
**Spain** +34 914 845 965

**Sweden/Norway/Finland**  
 +46 8 556 468 00  
**Switzerland** +41 61 48784 00  
**UK** +44 1442 233555  
**USA** +1 800 532 4752  
[www.thermo.com](http://www.thermo.com)

Thermo Electron Scientific Instruments LLC, Madison, WI USA is ISO Certified.  
 If the Class I option is not selected then the system is a Class IIIb laser product.



PS50445\_E 12/07M